


## Standard Form Crucibles



\*Reinforcing of the rim gives the crucible a more robust character enabling regular handling and reducing rim deformation.

Available Capacity: 2.5-700ml (larger sizes made to order).

Options: Lids, reinforced and/or bases.

|        |               |               |            | Estimated Weight (g) | Estimated Weight (g)  | Estimated Weight (g)   | Estimated Weight (g)    | Estimated Weight (g) |
|--------|---------------|---------------|------------|----------------------|-----------------------|------------------------|-------------------------|----------------------|
| Code   | Capacity (ml) | Diameter (mm) | Depth (mm) | Std. Form            | Rim Only (reinforced) | Base Only (reinforced) | Rim & Base (reinforced) | Std. Lid Only        |
| GCS2.5 | 2.5           | 12            | 15         | 3                    | 4                     | 4                      | 5                       | 2                    |
| GCS5   | 5             | 20            | 22         | 4                    | 5                     | 5                      | 6                       | 2                    |
| GCS8   | 8             | 25            | 24         | 6                    | 7                     | 7                      | 8                       | 2                    |
| GCS10  | 10            | 25            | 28         | 7                    | 7.5                   | 8                      | 9                       | 3                    |
| GCS15  | 15            | 30            | 31         | 11                   | 12                    | 13                     | 14                      | 4                    |
| GCS20  | 20            | 33            | 35         | 15                   | 16                    | 17                     | 18                      | 5                    |
| GCS25  | 25            | 36            | 34         | 19                   | 20                    | 21                     | 23                      | 5                    |

|        |     |     |     | Estimated Weight (g) | Estimated Weight (g) | Estimated Weight (g) | Estimated Weight (g) | Estimated Weight (g) |
|--------|-----|-----|-----|----------------------|----------------------|----------------------|----------------------|----------------------|
| GCS30  | 30  | 40  | 36  | 22                   | 24                   | 26                   | 27                   | 7                    |
| GCS35  | 35  | 41  | 38  | 25                   | 26                   | 27                   | 32                   | 7                    |
| GCS40  | 40  | 42  | 42  | 30                   | 32                   | 34                   | 36                   | 7                    |
| GCS50  | 50  | 45  | 44  | 38                   | 40                   | 42                   | 45                   | 8                    |
| GCS60  | 60  | 45  | 47  | 46                   | 48                   | 50                   | 54                   | 9                    |
| GCS70  | 70  | 49  | 53  | 50                   | 52                   | 54                   | 63                   | 12                   |
| GCS80  | 80  | 51  | 53  | 51                   | 54                   | 63                   | 72                   | 14                   |
| GCS90  | 90  | 54  | 57  | 52                   | 59                   | 71                   | 81                   | 15                   |
| GCS100 | 100 | 56  | 57  | 59                   | 68                   | 73                   | 90                   | 17                   |
| GCS110 | 110 | 56  | 60  | 68                   | 71                   | 81                   | 95                   | 19                   |
| GCS120 | 120 | 64  | 66  | 84                   | 88                   | 90                   | 98                   | 22                   |
| GCS150 | 150 | 65  | 69  | 100                  | 104                  | 107                  | 114                  | 21                   |
| GCS200 | 200 | 80  | 65  | 150                  | 159                  | 165                  | 180                  | 30                   |
| GCS250 | 250 | 81  | 73  | 175                  | 179                  | 183                  | 190                  | 32                   |
| GCS700 | 700 | 107 | 105 | 300                  | 305                  | 315                  | 320                  | 50                   |